

Medisoft 19 System Requirements

NOTE: Medisoft® 19 uses Advantage 11.10. If you are upgrading, you MUST install Advantage 11.10. Once you have installed Advantage 11.10, you will not be able to access an Advantage Server that is prior to Advantage 11. You will be able to access your database as long as your database server is at least Advantage 11.

Quick View: Medisoft 19 Supported Operating Systems

NOTE: Medisoft 19 is a 32-Bit application, and on a supported 64-Bit platform, the application will run in a 32-Bit mode.

OPERATING SYSTEM	SERVER		WORKSTATION		DATABASE
	Version of Medisoft 19		Version of Medisoft 19		Advantage 11.10
	Network Professional	Single-User or Advanced	Network Professional	Single-User or Advanced	Local or Remote
Windows Server 2003 SP2 32-Bit or 64-Bit, or Windows Server 2003 R2 SP2 32-Bit or 64-Bit	Y	N/A	Y	N/A	Remote
Windows Server 2012 or 2008 32-Bit	Y	N/A	Y	N/A	Remote
Windows Server 2012 or 2008 R2 64-Bit	Y	N/A	Y	N/A	Remote
Windows 8 Professional 32-Bit	Y	Y	Y	Y	Local
Windows 8 Professional 64-Bit	Y	Y	Y	Y	Local
Windows XP Professional SP3 32-Bit	N/A	Y	Y	Y	Local
Windows 8 or 7 Professional 32-Bit	Y	Y	Y	Y	Local
Windows 8 or 7 Professional 64-Bit	Y	Y	Y	Y	Local
Windows 7 Ultimate 32-Bit	Y	Y	Y	Y	Local
Windows 7 Ultimate 64-Bit	Y	Y	Y	Y	Local

Medisoft will not support Windows RT on the ARM processor.

Medisoft Hardware Requirements

Workstation (required)

CPU (Processor)	Intel Pentium 4 1.0GHz
RAM (Memory)	1GB
Storage Space	500MB (1.75GB for Net Pro)
Optical Drive	N/A
Network Card (NIC)	N/A
Display Monitor	1024x768 (1280x800 for widescreen displays)

Workstation (Recommended)

CPU (Processor)	Intel Core 2 Duo 1.6GHz
RAM (Memory)	2GB
Storage Space	1GB (1.75 for Net Pro)
Optical Drive	DVD-ROM
Network Card (NIC)	1Gbps

Server (Required)-Network Professional

CPU (Processor)	Intel Pentium 4 2.6GHz
RAM (Memory)	2GB
Storage Space	2GB
Optical Drive	DVD-ROM
Network Card (NIC)	1Gbps
Display Monitor	1024x768 (1280x800 for widescreen displays)

IMPORTANT: You also must have an Internet connection to use Network Professional.

Installation Overview

This document contains instructions for all versions of Medisoft 19 including Medisoft Reports Professional. Carefully review the instructions for your version of Medisoft 19 and hardware and software requirements before installing the software.

If you purchased Medisoft 19 Single User or Advanced, see [“Installation Instructions for Medisoft Single User and Advanced”](#) on page 5.

If you purchased Medisoft 19 Network Professional, see [“Installation Instructions for Medisoft Network Professional”](#) on page 9

If you purchased Medisoft Reports Professional, see [“Installation Instructions for Medisoft Reports Professional”](#) on page 17.

Only the installation of Medisoft Network Professional requires pre-installation setup of a network, as well as installation on both the server computer and the workstations.

If you receive an error trying to open a PDF file or to install Adobe Reader, go to the Adobe web site and download the latest version of Adobe Reader.